

Grid Box

#54 904

T

19. Glutaraldehyde (3%) - 5 min. then CCl_4 /Osmium(80%) - 2 min. then
Glutaraldehyde (3%) -10 min.

U

20. Glutaraldehyde (3%) - 15 min. then CCl_4 / Osmium (80%) - 5 min.

Grid Box

#54 905

A

21. CCl_4 / Osmium (80%) - 2 min. then Glutaraldehyde (3%) - 15 min.

B

22. CCl_4 /Osmium (80%) - 2 min. then Glutaraldehyde (3%) - 10 min. then
Osmium (2%) - 10 min.

23. CCl_4

C

23. CCl_4 / Osmium (80%) - 2 min.

D

24. Glutaraldehyde (3%) - 5 min. then CCl_4 / Osmium (80%) - 2 min.

